

The following **16** technologies across **eight** newborn product categories are currently qualified by NEST360°.

Product Category	Product Name*
Phototherapy	MTTS Colibri Phoenix Brilliance Pro
Glucometer	Nova StatStrip Glucose Hospital Meter System Nova StatStrip Xpress2
CPAP	Diamedica UK Baby CPAP 10 Pumani bubbleCPAP
Flow Splitter	Caire Chart SureFlow Oxygen Flow Station Canta Flow Splitter
Oxygen Concentrator	Canta V8-WN-NS Caire Airsep Newlife Intensity 10 Longfian DoubleFlow
Pulse Oximeter	Acare Lifebox Bistos BT-710
Suction Pump	3A Aspeed Professional
Radiant Warmer	MTTS Wallaby Phoenix NWS-101

*Product names listed in alphabetical order

What is “NEST360° qualified” and how are technologies qualified by NEST360°?

NEST360° evaluates technologies for newborn care in low-resource settings using the five step process described below. Technologies are considered NEST360° qualified if they meet pre-defined performance metrics in each step.

Step 1 – Identification:

NEST360° performs desk research to identify newborn technologies that are currently on the market or that are in development across 16 product categories^[1]. This includes reviewing technologies listed in the [Newborn Technology Landscape](#), which is updated every six months.

Step 2 – Alignment with Target Product Profile:

Technologies are compared against the [Newborn Target Product Profiles \(TPPs\)](#), matching characteristics defined in the TPPs to the comparable characteristics outlined in the manufacturer’s package inserts and other publicly available materials. Each technology is ranked according to the number of characteristics that meet optimal or minimal TPP characteristics. Lead candidates within each product category are based on the rankings. To date, an initial set of eight product categories^[2] have multiple commercially available options that meet a majority of TPP characteristics.

Step 3 – Technical Testing:

Technologies that meet a majority of TPP characteristics are then purchased for laboratory measurement of each characteristic in the TPP. To date, two or three technologies were measured within each of the eight product categories^[2], with two units of each technology measured by three unique users.

Step 4 – Environmental Testing:

Environmental testing is performed on the qualifying technologies exposing the technologies to extreme heat, humidity, and dust levels as well as voltage surges and sags in order to determine which technologies perform under these conditions. Functionality test protocols, which compare performance to baseline technical measurements, are used to assess technologies at discrete points in time throughout each environmental test.

Step 5 – Usability Testing:

In parallel with technical and environmental testing, comparative usability testing is conducted on each technology with at least five unique users (nurses) in the United States and in Malawi, respectively. Users are asked to perform a series of tasks with each technology after watching a brief instructional video. Metrics measured include efficiency, effectiveness, and satisfaction, with each user ultimately identifying a preferred technology.

^[1] 16 Product Categories: Syringe Pump, Bilirubinometer, Phototherapy, Glucometer, Hemoglobinometer, pH Monitor, Sepsis Diagnostic, CPAP, Flow Splitter, Oxygen Concentrator, Pulse Oximeter, Respiratory Rate / Apnea Monitor, Suction Pump, Radiant Warmer, Temperature Monitor, Conductive Warmer

^[2] Eight Product Categories with devices that are currently NEST360° qualified: Phototherapy, Glucometer, CPAP, Flow Splitter, Oxygen Concentrator, Pulse Oximeter, Suction Pump, Radiant Warmer